



Letter of Declaration

To whom it may concern,

We, HANGZHOU ALLTEST BIOTECH CO., LTD., located on #550, Yin Hai Street, Hangzhou Economic & Technological Development Area, Hangzhou -310018, P.R. China

Declare that:

SARS-CoV-2 Antigen Rapid Test (COVID-19 Antigen Rapid Test) is designated to detect nucleocapsid protein (NP) of SARS-CoV-2 virus. Epitope of capture antibody and detection antibody was not in the mutation area. Thus, above mentioned spike protein mutation SARS-CoV-2 and nucleocapsid protein mutation SARS-CoV-2 of Alpha(B.1.1.7), Beta(B.1.351), Gamma(P.1), Delta(B.1.617.2), Lambda(C.37), Zeta(P.2), Kappa(B.1.617.1), Eta(B.1.525), Iota(B.1.526), Epsilon(B.1.427), Theta(P.3) and Omicron(B.1.1.529) is not affecting virus performance of SARS-CoV-2 Antigen Rapid Test (COVID-19 Antigen Rapid Test).

Signature:

R&D Manager: Kael Chen

Date: Nov. 28th, 2021

Hangzhou Alltest Biotech CO., Ltd

Add:#550 Yin Hai Street, Hangzhou Economic and Technologic Development Area
310018, P.R China

Оригинальный перевод:

Всем заинтересованным лицам

На запросы относительно мутации вируса SARS-CoV-2 (COVID-19) специфичности и чувствительности зарегистрированных экспресс-тестов.

Мы, Ханчжоу Аллтест Биотек КО., ЛТД., Расположены на #550, Улица Инхай, Район экономического и технологического развития Ханчжоу, Ханчжоу -310018, П.Р. Китай

Заявляем, что:

Экспресс Тест на Антиген SARS-CoV-2 (COVID-19 Antigen Rapid Test) предназначен для выявления нуклеокапсидного белка (НБ) SARS-CoV-2 вируса. Эпитоп антитела захвата и антитела обнаружения не находился в области мутации. Таким образом, вышеуказанная мутация белка шипа SARS-CoV-2 и мутация белка нуклеокапсида SARS-CoV-2 из Альфа(В.1.1.7), Бета(В.1.351), Гамма(Р.1), Дельта(В.1.617.2), Лямбда(С.37), Зета(Р.2), Каппа(В.1.617.1), Эта(В.1.525), Йота(В.1.526), Эпсилон(В.1.427), Тета(Р.3) и Омикрон(В.1.1.529) не влияет на вирусную активность SARS-CoV-2 Экспресс Теста на Антиген (COVID-19 Antigen Rapid Test).

Управляющий НИОКР: Kael Chen

Дата: Ноябрь, 2021 Ханчжоу АллТест Биотек

Приложение:#550 Улица Инхай, Район экономического и технологического развития Ханчжоу
310018, Китайская Народная Республика